

SelectaDNA High Velocity Tagging System

Irrefutably identifying an individual in a crowd or at distance has always presented law enforcement officials with a number of challenges. The SelectaDNA High Velocity Tagging System allows a DNA coded pellet to be used to mark an individual, from a safe distance, so that the individual can be identified as being present at a specific location at a specific time.

- On contact with the target the uniquely coded SelectaDNA solution will leave a unique DNA code on the target to provide evidence that will enable the relevant authorities to confirm or exclude a person or object from their involvement in a situation under scrutiny.
- The SelectaDNA High Velocity delivery system will enable Police, Military and other Government Law Enforcement Agencies personnel to remain at a safe distance (up to approximately 30-40 metres) from persons of interest while deploying the SelectaDNA coded pellet.
- The SelectaDNA High Velocity delivery system comes in two formats; either a kit which allows the launcher to be configured as a rifle or a pistol to suit the requirements of the intended application or alternatively as a pistol.
- The delivery systems are powered by a 12g powerlet offering up to 20 shots per powerlet. The powerlets are single use items and must be used on the day of first activation/ connection and removed from the delivery system at the end of a deployment.
- The SelectaDNA pellets are loaded into 8 shot magazines for deployment.
- The SelectaDNA pellets are supplied in packs of 16 pellets. All pellets in a pack are filled with the same uniquely coded DNA and UV solution to allow the location of the DNA code to be identified under UV light.



For more information visit

www.selectadna.co.uk

 facebook.com/SelectaDNA

or contact

+44 (0)1689 860757

Selectamark Security Systems plc, 1 Locks Court,
429 Crofton Road, Locksbottom, Kent, BR6 8NL. UK.

 @SelectaDNA #SelectaDNA



SelectaDNA is a registered trademark of:



SelectaDNA High Velocity Tagging System



Description	PISTOL	RIFLE
Product Code:	19410	19420
Description:	Plastic body and anodized barrel. Semi automatic, open bolt blow forward, 1 trigger pull = 1 shot (6 per second)	Plastic body and anodized barrel. Semi automatic, open bolt blow forward, 1 trigger pull = 1 shot (6 per second)
Dimensions:	280mm x 204mm high/11 Inches x 8 Inches high	840mm x 230mm/33 Inches x 9 Inches high
Calibre	.68	.68
Colour:	Orange/black	Orange/black
Unloaded weight :	1.23kg/2.7lbs	2.2kg/4.85lbs
Pellet Capacity:	8 shot magazine	8 shot magazine
SelectaDNA Pellets (Part Number 19430):	16 single code DNA Pellets	16 single code DNA Pellets
Power Source:	12g powerlet	12g powerlet
Power source capacity*	20 shots	20 shots

* Power sources must be used on day of first activation/connection
Product specifications are subject to change at the manufacturers discretion.

For more information visit

www.selectadna.co.uk

 facebook.com/SelectaDNA

or contact

+44 (0)1689 860757

Selectamark Security Systems plc, 1 Locks Court,
429 Crofton Road, Locksbottom, Kent, BR6 8NL. UK.

 @SelectaDNA #SelectaDNA



SelectaDNA is a registered trademark of:



Selectamark Secure Asset Register
Certificate 0123-10003
Certificate Number 012345

Assessed to ISO 9001:2008
Certificate Number 012345

IS 602150